



Prontosan®
Askina® Calgitrol®

WOUND BED PREPARATION TAKEN SERIOUS

REMOVE BARRIERS TO HEALING

BIOFILMS & WOUND INFECTIONS

A BIOFILM IS A COMPLEX STRUCTURE OF BACTERIA EMBEDDED IN A PROTECTIVE SLIMY MATRIX OF SUGARS AND PROTEINS.

HOW BIOFILMS DELAY WOUND HEALING¹⁷

- Induce a chronic inflammatory response
- Increase level of exudate which in turn perpetuates the biofilm
- Protect the bacteria from antibiotics, antimicrobials as well as host defenses
- Impair migration and proliferation of keratinocytes

ACUTE AND CHRONIC WOUNDS HAVE A TENDENCY TO GET INFECTED AND COMPROMISE THE NORMAL HEALING PATHWAY, LEADING TO A GREATER BURDEN ON HEALTH SYSTEMS, LONG TERM DISABILITIES AND AN OVERALL REDUCTION OF A PATIENT'S QUALITY OF LIFE.



BIOFILM

90 % of chronic wounds contain biofilm with a role in wound infection.¹



INFECTION

50 % of chronic wounds infected²
Half of chronic wounds are estimated to be infected.



PAIN

50 % of chronic wound patients have painful wounds
Pain is present in almost half of patients with chronic wounds.³



NURSING TIME

Up to 50 % of nursing time⁴
Almost half of nursing time is dedicated to the management of chronic wounds.



IMPACT

3 % of total health expenditures⁵
Chronic wounds cause pain, loss of function and mobility, depression, distress and anxiety, embarrassment, social isolation and prolonged morbidity.



ESTIMATED COST

Up to 10 times more cost for complications⁶
Cost of healing increase as the time to heal is prolonged and the incidence of complications is higher.

CLEAN HANDS SAVE LIVES

HAND HYGIENE IS THE FIRST MEASURE TO PREVENT INFECTION

HAND HYGIENE GUIDELINES

When to wash?

- When hands are visibly dirty or soiled with blood or other body fluids.
- After using the toilet.
- If exposure to spore forming pathogens is suspected or confirmed, i.e. Clostridium Difficile.

All other instances, use ABHR for routine decontamination of hands.

50%

of all health care-associated infections (HAI) can be avoided by proper hand hygiene.¹²

WHY SELECT AN ALCOHOL BASED HAND RUB?

- Can be used everywhere.
- Faster than hand washing.
- More effective than washing.
- More gentle to your skin.


TAKE NOTE:

- Reduce the bacterial load on the patient's skin by decolonization.
- Disinfect surfaces and instruments.

5 MOMENTS FOR HAND HYGIENE IN WOUND CARE

IMMEDIATELY BEFORE: Touching the wound dressing / site

- Before donning gloves
- Before physically examining the wound site, including before taking wound swabs.
- Before removing stiches / clips from post-operative wounds.
- Before preparing the necessary items for changing the wound dressing.
- Before changing the wound dressing.



IMMEDIATELY AFTER: Any task involving potential body fluid exposure

- After removing gloves
- After wound examination / sample collection.
- After removing stiches / clips from post-operative wounds.
- After changing a wound dressing.



ANTISEPTIC SCRUB

Bioscrub

4 % Chlorhexidine Gluconate



ALCOHOL-BASED HAND RUBS

D-Germ

70 % Propyl Alc, 0.5 % CHG and emollient (Passed EN1500/12791/13624/13727)¹⁶

D-Bac

75 % Isopropyl with emollient WHO II Formulation

Softa-Rub

Ethanol 78.1 %, Isopropyl 7.81 % and emollients

Softa-Gel

Ethanol 78.1 %, Isopropyl 3.8 % and emollients

Prontosan® Wound Irrigation Solution

Prontosan® Wound Gel I Wound Gel X

WOUND BED PREPARATION TAKEN SERIOUSLY

INDICATIONS

Prontosan® Wound Irrigation Solution is indicated for cleansing and moistening of acute and chronic wounds.

Prontosan® Gel and Gel X act as an effective barrier to reduce microbial penetration through the dressing and to decontaminate the wound bed between dressing changes.

MODE OF ACTION

Unique combination of Betaine: A gentle effective surfactant to penetrate, clean and remove wound debris and biofilm **and Polihexanide:** Acts as preservative and reduces the risk of infections.

ADVANTAGES

- Management and prevention of biofilm reformation.^{7,8}
- Helps to prevent infections.⁹
- Improved patient outcomes, including time to heal.¹⁰
- Well-known substances with low allergenic potential.¹¹
- Can be used up to 8 weeks after first opening.
- Compatible with commonly used wound care products.
- Complimentary to Acelyty VAC Ultra.



Prontosan® Wound Irrigation Solution may be used to saturate compresses, gauzes or pads for cleansing as needed. Prontosan should be used with every dressing change in order to achieve and maintain a visually clean wound.



For application in deep or tunneling wounds, wound cavities and difficult to access areas, apply a 3–5 mm layer of Prontosan® Wound Gel and cover with a secondary dressing.



For superficial or larger surface area wounds, apply a 3–4 mm layer of the highly viscous Prontosan® Wound Gel X and cover with a secondary dressing.

Prontosan® Debridement Pad

SOFT MECHANICAL REMOVAL OF SLOUGH AND DEBRIS


INDICATIONS

Prontosan® Debridement Pad has been designed to support the Wound Bed Preparation with Prontosan® Wound Irrigation Solution. Prontosan® Debridement Pad frees the wound from coatings and dead cell residues (debris) and absorbs excess exudate and slough. Intact tissue is spared.

Prontosan® Debridement Pad produces good results even with scaly and necrotic coatings, if they are subject to prior autolytic treatment.


HOW DO MICROFIBERS WORK

Microfibers are able to attach themselves to even the smallest, most microscopic particles, ones that normal monofilament fibers (positively giant in comparison) crudely brush past. If adhesive forces were visible, you would be able to see them between the microfibers and particles. There are millions of microfibers in a pad, so the overall sticking effect is magnified dramatically. That's why cell residues (debris), exudate and slough can be „hoovered“ up by microfibers.¹⁵



APPLICATION

The unmarked side of the pad has to be moistened with 15 – 20 ml of Prontosan® Wound Irrigation Solution. Apply gentle pressure, wipe the moistened side over the wound and / or adjacent skin in a circular motion.



Following debridement, it is advisable to irrigate the wound thoroughly for a second time with Prontosan® Wound Irrigation Solution and to continue with standard treatment.

ADVANTAGES

- Good cleansing and debridement due to microfiber technology.
- Rapid and soft debridement, no tissue irritation.
- Unique droplet shape to allow debridement of cavities and areas difficult to reach.
- Blister packaging to allow safe and aseptic soaking of the pad prior to use.
- Flexible, thin debridement sheet consisting of a cleansing layer (Polyester / Polypropylene) and an absorbent backing layer (Polypropylene).



PRONTOSAN® DEBRIDEMENT PAD REDUCE TIME

When combining Prontosan Solution® with the pad, soaking time is reduced.

Askina® Calgitrol® Range

PATENTED IONIC SILVER ALGINATE MATRIX

INDICATIONS

Infected or critically colonized wounds:

- Pressure ulcers
- Venous and arterial ulcers
- Diabetic foot ulcers
- Second degree burns
- Traumatic wounds
- Donor sites

ADVANTAGES

- 100 % ionic silver, active against a broad spectrum of microorganisms.¹³
- Immediate availability of silver ions.¹⁴
- Sustained controlled release of silver ions during use of the dressing.¹⁴
- Non staining.
- Different presentations adapted to the wound bed based on the same technology.



Askina® Calgitrol® Ag

A sterile dressing consisting of two layers:

- An absorbent polyurethane foam layer for the absorption of wound exudate.
- An ionic silver alginate matrix for broad antimicrobial effectiveness and helps reduce the bacterial load.



Askina® Calgitrol® THIN


A thin layer of ionic silver alginate matrix. Soft and conformable, it is well adapted for deep and difficult-to-dress wounds.



Askina® Calgitrol® Paste

Its conformability allows an extremely close contact between the active ionic silver alginate matrix and the wound bed, which is particularly valuable in difficult to manage wounds such as tunnels and sinuses.

MODE OF ACTION




Silver dressings using Calgitrol® technology are ready for use and do not require prior wetting to prime activity.

Moisture absorption through the matrix leads to softening and swelling of the alginate structure, which facilitates controlled and steady state release of ionic silver.

B.Braun Wound Care Products

TECHNOLOGIES TO SUPPORT WOUND HEALING AT EVERY PHASE

NECROTIC	DRY FIBRIN	MOIST FIBRIN	GRANULATION	EPITHELIALIZATION	INFECTION
					
Askina® Gel Hydrogel		Askina® Sorb Alginate		Askina® Biofilm Transparent Hydrocolloid	Askina® Carbosorb Charcoal dressing
					
		Askina® Silnet Silicone Wound Contact Layer			Biocream PVP-I Cream
					
		Askina® Foam Foam			Askina® Calgitrol AG Silver Alginate Foam
					
			Askina® Dressil Silicone Foam		Askina® Calgitrol Thin Silver Alginate Matrix
					
			Askina® Dressil Border Silicone Foam with Border		Askina® Calgitrol Paste Silver Alginate Paste
					
		Prontosan® Irrigation Solution Betaine Et Polihexanide Wound Cleanser			
		Prontosan® Debridement Pad** Microfiber Debridement Pad			
		Prontosan® Wound Gel / Gel X Betaine Et Polihexanide Wound Gel			

**Can be used on necrotic wounds if prior autolytic treatment is implemented

Product Code	Product Description	Size	Nappi Code
WOUND CLEANSERS			
400484	PRONTOSAN SOLUTION	40 ML	714724-001
400403	PRONTOSAN SOLUTION	350 ML	716072-001
400446	PRONTOSAN SOLUTION	1 L	718761-001
400240A	PRONTOSAN SOLUTION + CAPS/ADAPTOR	10 x (1 L + CAPS)	720609-001
400565	PRONTOSAN WOUND SPRAY	75 ML	721428-001
400505	PRONTOSAN WOUND-GEL	30 ML	715025-001
400517	PRONTOSAN WOUND-GEL X	50 GR	718945-001
400508	PRONTOSAN WOUND-GEL X	250 GR	718945-002
3908456	PRONTOSAN DEBRIDEMENT PAD	3 / PACK	288951-001
3908457	PRONTOSAN DEBRIDEMENT PAD	10 / PACK	288951-002

Product Code	Product Description	Size	Nappi Code
THIN HYDROCOLLOID			
F72090	ASKINA BIOFILM TRANSPARENT	10X10 CM	415959-001
F72097	ASKINA BIOFILM TRANSPARENT	15X15 CM	415965-001
F72098	ASKINA BIOFILM TRANSPARENT	20X20 CM	415966-001

Product Code	Product Description	Size	Nappi Code
HYDROGEL			
001419N	ASKINA GEL	15 GR	415957-001
14291	ASKINA GEL	100 GR	230325-001

Product Code	Product Description	Size	Nappi Code
SILICONE WOUND CONTACT LAYER			
5197510	ASKINA SILNET	10X7.5 CM	113650-001
5191810	ASKINA SILNET	10X18 CM	113648-001
5192305	ASKINA SILNET	20X30 CM	113649-001

Product Code	Product Description	Size	Nappi Code
SILICONE FOAM			
5295710	ASKINA DRESSIL	5X7 CM	204195-001
5291010	ASKINA DRESSIL	10X10 CM	204196-001
5291210	ASKINA DRESSIL	10X20 CM	204197-001
5291510	ASKINA DRESSIL	15X15 CM	204198-001
5292005	ASKINA DRESSIL	20X20 CM	204199-001

Product Code	Product Description	Size	Nappi Code
SILICONE FOAM WITH BORDER			
5396610	ASKINA DRESSIL BORDER	6X6 CM	204179-001
5397510	ASKINA DRESSIL BORDER	7.5X7.5 CM	204180-001
5391010	ASKINA DRESSIL BORDER	10X10 CM	204181-001
5391210	ASKINA DRESSIL BORDER	10X20 CM	230324-001
5391510	ASKINA DRESSIL BORDER	15X15 CM	204182-001
5392005	ASKINA DRESSIL BORDER	20X20 CM	204185-001
5395210	ASKINA DRESSIL BORDER	15X20 CM	204183-001
5491605	ASKINA DRESSIL SACRUM	16X17.5 CM	204184-001
5492110	ASKINA DRESSIL SACRUM	21X22 CM	230322-002
5592205	ASKINA DRESSIL HEEL	22X21.6 CM	230323-001

Product Code	Product Description	Size	Nappi Code
FOAM			
7240105	ASKINA HEEL	225 CM²	526956-001
7240710	ASKINA FOAM	5X7 CM	526958-001
7241010	ASKINA FOAM	10X10 CM	440942-001
7241210	ASKINA FOAM	10X20 CM	440943-001
7242005	ASKINA FOAM	20X20 CM	526964-001
7244010	ASKINA FOAM CAVITY	2.5X40 CM	526966-001

Product Code	Product Description	Size	Nappi Code
ALGINATE			
2100S	ASKINA SORB	6X6 CM	527463-002
2101S	ASKINA SORB	10X10 CM	527462-002
2102S	ASKINA SORB	15X15 CM	527465-001
2105S	ASKINA SORB ROPE	2.7X34 CM	527469-001

Product Code	Product Description	Size	Nappi Code
SILVER ALGINATE			
6201010	ASKINA CALGITROL THIN	10X10 CM	113646-001
6202010	ASKINA CALGITROL THIN	20X20 CM	113647-001
6202403	ASKINA CALGITROL THIN	20X40 CM	139998-002
6211010	ASKINA CALGITROL AG	10X10 CM	440944-001
6211510	ASKINA CALGITROL AG	15X15 CM	440946-001
6212010	ASKINA CALGITROL AG	20X20 CM	440948-001
6241505	ASKINA CALGITROL PASTE	15 GR	169686-001
6241001	ASKINA CALGITROL PASTE	100 GR	230321-001
6242501	ASKINA CALGITROL PASTE	250 GR	274233-001
6245001	ASKINA CALGITROL PASTE	500 GR	274234-001

Product Code	Product Description	Size	Nappi Code
ODOUR MANAGEMENT			
9025006	ASKINA CARBOSORB	10X10 CM	646529-001
9025014	ASKINA CARBOSORB	10X20 CM	646530-001

Product Code	Product Description	Size	Nappi Code
PVP-IODINE CREAM			
DM0501398	BIOCREAM	25 GR	871885-018
DM0501398T	BIOCREAM TENDER	25 GR	-
DM0501388	BIOCREAM	500 GR	871885-026

Product Code	Product Description	Size	Nappi Code
HAND CREAM			
18605	TRIXOLIND	100 ML	-

Product Code	Product Description	Size	Nappi Code
DRESSING TRAYS			
BBMZA70101	DRESSING TRAY	SMALL	201573-001
BBMZA70102	DRESSING TRAY	LARGE	201574-001

Product Code	Product Description	Size	Nappi Code
ANTISEPTIC SCRUB			
DM0501296	BIOSCRUB	50 ML	871893-010
DM0501297	BIOSCRUB	500 ML	871893-002
DM0501297P	BIOSCRUB WITH PUMP	500 ML	715648-001

Product Code	Product Description	Size	Nappi Code
ALCOHOL-BASED HAND RUBS/GELS			
DM0501308P	D-GERM WITH PUMP	500 ML	871907-003
DM0501464	D-BAC WITH PUMP	500 ML	555522-003
DM0501472P	SOFTA-RUB WITH PUMP	500 ML	555522-003
DM0501470P	SOFTA-GEL WITH PUMP	500 ML	555522-003

Product Code	Product Description	Size	Nappi Code
SURFACE DISINFECTANT			
DM0501100	ALCOHOL SURFACE DISINFECTANT	500 ML	555522-003
19894	MELISEPTOL WIPES SENSITIVE	100/ FLOW PACK	-

Product Code	Product Description	Size	Nappi Code
HIGH-LEVEL DISINFECTANT (INSTRUMENTS/ ENDOSCOPES)			
19812	STABIMED ULTRA	800 G	-
19939	STABIMED ULTRA	4K G	-

B. Braun Medical (Pty) Limited | P.O. Box 1787, Randburg, 2125 | South Africa
 Tel. 010 222 3000 | Fax 086 230 0517 / 010 222 3066 | Registration No.1995/003058/07 | www.bbraun.com

REFERENCES

- Clinically addressing biofilm in chronic wounds: Attinger C, Wolcott R, *Advances in Wound Care* 2012; 1(3):127-132
- C. Dowsett: Adopting the two-week challenge in practice: making the case for silver dressings, *Wounds UK*, Vol. 10, N°2, 2014
- J Guest et al. Health economic burden that wounds impose on the National Health Service in the UK. *BMJ Open*, 2015;5:e009283. doi: 10.1136/bmjopen-2015-009283
- H Edwards et al. Health service pathways for patients with chronic leg ulcers: identifying effective pathways for facilitation of evidence based wound care. *BMC Health Services Research*, 2013; 13(86), <https://doi.org/10.1186/1472-6963-13-86>
- K Järbrink et al. Prevalence and incidence of chronic wounds and related complications: a protocol for a systematic review. *Systematic Reviews*, 2016; 5(1):152
- G. Bennett, C. Dealey and J. Posnett: The cost of pressure ulcers in the UK, *Age and Aging*, 33, 2004
- Efficacy of various wound irrigation solutions against biofilms. Seipp HM, Hofmann S, Hack A, Skowronsky A, Hauri A., *ZfW* 2005;4(5):160-163.
- Davis SC, Harding A, Gil J, Parajón F, Valdes J, Solis M & Higa A "Effectiveness of a polyhexanide irrigation solution on methicillin-resistant *Staphylococcus aureus* biofilms in a porcine wound model" in *International Wound Journal* ISSN 1742-4801, 2017, 1-8. © 2017 Medicalhelplines.com Inc and John Wiley & Sons Ltd doi: 10.1111/iwj.12734.
- Moore, M 0.1% Polyhexanide-Betaine Solution as an Adjuvant in a Case-Series of Chronic Wounds, *Surg Technology International*, 2016.
- Bellingeri, A. et al. "Effect Of A Wound Cleansing Solution On Wound Bed Preparation And Inflammation In Chronic Wounds: A Single-Blind RCT". *Journal of Wound Care* 25.3 (2016): 160-168. Web.
- Evaluation of the efficacy and tolerability of a solution containing propyl betaine and polyhexanide. Romanelli M, Dini V, Barbarera S, Bertone MS. *Skin Pharmacol Physiol* 2010;23(Suppl 1):41-44.
- Pittet, Didier; „Adapt to adopt; TEDxPlaceDesNations“, under <https://www.youtube.com/watch?v=5tgH0uIqqc> (accessed at 2 May 2016)
- Instruction for use: Askina® Calgitrol® Ag, Askina® Calgitrol® THIN, Askina® Calgitrol® Paste.
- Opasanon S, Magnette A, Meuleniére F, Harding K. Askina® Calgitrol® Made Easy. *Wounds International* 2012; 3(1). Available from www.woundsinternational.com
- Prontosan Debridement Pad Made Easy, *Wounds International*, May 2018
- EN1500/EN12791/EN13624/EN13727 January 2018
- Biofilm made easy. Vol. 1 Issue 3, May 2010 (www.woundsinternational.com/pdf/content_8851.pdf)

For full prescribing information, please refer to the package insert approved by the medicine regulatory authority.