




AIS CONSUMABLES


Auto Infusion Systems:

- Infusomat SpacePlus Dedicated Sets
- Infusomat Space Dedicated Sets
- Infusomat CompactPlus Dedicated Sets
- Perfusor SpacePlus Syringes
- Perfusor Space Syringes
- Perfusor CompactPlus Syringes
- IV Accessories

Infusomat® Space Dedicated Sets

	Product Code	Nappi Code	Description	Unit of Sale
	8700036SP	437403-001	INFUSOMAT SPACE LINE STANDARD, NO Y-Site, 250cm	EACH
	8270350SP	150353-001	INFUSOMAT SPACE LINE STANDARD, NO Y-Site, 300cm	EACH
	8700087SP	484118-003	INFUSOMAT SPACE INJECTION TUBING with Injection Port, PVC, 250cm	EACH

FILTERED SETS

	8250731SP-TPN	160255-001	INFUSOMAT SPACE LINE NEUTRAPUR PVC-FREE with 1.2 Micron Fliter (incl TPN), NO Y-SITE, 250cm	EACH
	8700095SP	156401-001	INFUSOMAT SPACE LINE NEUTRAPUR, PVC FREE with 0.2 micron filter, no y-site, 250cm	EACH

BLOOD SETS

	8270066SP-26	643688-005	INFUSOMAT SPACE LINE FOR TRANSFUSION WITH BLOOD FILTER (200 µm pore size), 250cm	EACH
---	--------------	------------	--	------

UV PROTECT SETS

	8700128SP	210130-001	INFUSOMAT SPACE LINE, AMBER UV-PROTECTED, 250cm, Neutrapur	EACH
---	-----------	------------	--	------

LOW ABSORPTION SETS

	8250731SP	693001-001	INFUSOMAT SPACE LINE, NEUTRAPUR PVC-FREE, ADMIN OF CYTOTOXIC DRUGS, 250cm	EACH
---	-----------	------------	---	------

Infusomat® Space Dedicated Sets

	Product Code	Nappi Code	Description	Unit of Sale
	8270350SPNF-NICU	160249-001	INFUSOMAT SPACE LINE with 2 Needle Free Valves and 0.2 Micron Filter, 300cm	EACH
	8700110SP	604052-001	INFUSOMAT SPACE LINE with 1 Needle Free Valve, PUR TUBE, 250cm	EACH
	8270350SPNF	461675-001	INFUSOMAT SPACE LINE with 2 Needle Free Valves, PVC, 300cm	EACH
	8270350SPS	461674-001	INFUSOMAT SPACE LINE with 1 Needle Free Valve, 300CM	EACH
	8270055SPNF	447904-001	INFUSOMAT SPACE LINE with 2 Needle Free Valves, 250CM	EACH
	8700036SPS	447902-001	INFUSOMAT SPACE LINE with 1 Needle Free Valve, 250CM	EACH

FEEDING SETS

	8250857SP	597529-001	INFUSOMAT SPACE LINE, ENTERAL NUTRITION, (purple connectors), 250cm	EACH
	8250832SP	1052317-001	INFUSOMAT SPACE LINE, ENFit, PVC Free, 250cm	EACH

EPIDURAL SETS


	8700131SP	1001985-001	INFUSOMAT SPACE LINE, Epidural, PVC Free, 250cm	EACH
---	-----------	-------------	---	------

BURETTES (DOSIFIX)

	8700141SP	1120043-001	INFUSOMAT SPACELINE, needle free valve, DOSIFIX, PVC Free, 270 cm	EACH
	8700140SP	1120044-001	INFUSOMAT SPACELINE DOSIFIX, PVC Free, 270 cm	

Infusomat® Plus Dedicated Sets

	Product Code	Nappi Code	Description	Unit of Sale
	8700210	1048145-001	INFUSOMATPLUS LINE, drip chamber incl. AirStop function, 15 micron fluid filter	EACH
	8700250	1048146-001	INFUSOMATPLUS LINE, needle free, drip chamber incl. AirStop function, 15 micron fluid filter, PVC free	EACH
	8700360	1048147-001	TRANSFUSION SET, 200 micron blood filter, DEHP free, spinlock patient connector	EACH
	8700270	1048148-001	UV LIGHT PROTECTED SET, 15 micron fluid filter, drip chamber with airstop function, luer lock patient connector	EACH
	8700280	1048149-001	FLUSHING SET, 15 micron fluid filter, inline drip chamber with airstop function, spinlock patient connector	EACH



Perfusor® Space Syringes

	Product Code	Nappi Code	Description	Unit of Sale
	8728810F-06	502065-001	PERFUSOR SYRINGE, 50ML CLEAR, WITH NEEDLE, Luer Lock	EACH
	8728861F-06	210134-001	PERFUSOR SYRINGE, 50ML AMBER, WITH NEEDLE, Luer Lock	EACH
	4617509F	677493-003	OMNIFIX SYRINGE, 50ML, CLEAR WITHOUT NEEDLE, Luer Lock	EACH
	4616026-01	1052312-001	ENTERAL FEEDING SYRINGE 5ml, purple, ENFit, PVC free	EACH
	4616027-01	1052313-001	ENTERAL FEEDING SYRINGE 10ml, purple, ENFit, PVC free	EACH
	4616028-01	1052314-001	ENTERAL FEEDING SYRINGE 20ml, purple, ENFit, PVC free	EACH
	4616029-01	1052315-001	ENTERAL FEEDING SYRINE 50ml, purple, ENFit, PVC free	EACH

SYRINGE PACKS

	BBMZA8728828F-OPCK	168436-001	PERFUSOR PACK, 50ml Amber Omnifix Syringe & 150cm Ext Set	EACH
	BBMZA8728844F-CPCK	168437-001	PERFUSOR PACK, 50ml Clear Omnifix Syringe & 150cm Ext Set	EACH
	BBMZA4617029V	199555-001	PERFUSOR PACK, 2ML Omnifix Syringe & 150cm Ext	EACH
	BBMZA4617053V	199556-001	PERFUSOR PACK, 5ML Omnifix Syringe & 150cm Ext	EACH
	BBMZA4617100V	199557-001	PERFUSOR PACK, 10ML Omnifix Syringe & 150cm Ext	EACH
	BBMZA4617207V	199558-001	PERFUSOR PACK, 20ML Omnifix Syringe & 150cm Ext	EACH

	Product Code	Nappi Code	Description	Unit of Sale
EXTENSION SETS				
	8722960	462098-001	PERFUSOR TUBING N 150CM	EACH or BOX/50
	8722935	462128-008	PERFUSOR TUBING P.E. 150CM	EACH or BOX/50
	8723017	210424-001	PERFUSOR TUBING AMBER, LIGHT SENSITIVE 150CM	EACH or BOX/50
	8726019	505960-001	PERFUSOR TUBING, PCA - WITH BACK-CHECK VALVE,	EACH or BOX/50
	471950	452099-004	T-PORT EXTENSION SET ET-04T	EACH or BOX/50
	4092945	551379-001	SPIRAL TUBING - 300CM	EACH or BOX/50
	87229910	1052316-001	ORIGINAL PERFUSOR ENFit EXTENSION SET	EACH or BOX/50



	Product Code	Nappi Code	Description	Unit of Sale
CONES, STOPPERS, ADAPTORS				
	4495101	419753-001	COMBI-RED MALE/FEMALE CLOSING LUER LOCK	EACH or BOX/100
	4097076	418773-009	CLOSING CONES, LUER LOCK	EACH or BOX/100
	4238010	446882-003	IN-STOPPERS, INJECTION PORT, LATEX FREE	EACH or BOX/100
	415020	484575-003	FILTER STRAW	EACH or BOX/100
	5206642	555522-003	ADAPTOR MALE/MALE - CONNECTORS	EACH or BOX/100
	4095111	633950-007	DISCOFIX 3 WAY STOPCOCK BLUE ■	EACH or BOX/100
	4095120	633951-007	DISCOFIX 3 WAY STOPCOCK RED ■	
	4095146	633952-007	DISCOFIX 3 WAY STOPCOCK WHITE □	
	4099753	634175-006	INTRAPUR PAED I.V. FILTER 0.2 MICRON	EACH or BOX/100
	4099850	521572-006	INTRAPUR PAED LIPID I.V.FILTER 1.2 MICRON	
	4099702	634173-006	INTRAPUR LIPID I.V. FILTER 1.2 MICRON	
	4099800	634176-006	INTRAPUR PLUS I.V. FILTER 0.2 MICRON	



